



Setup details

Unistat® 405w & Glas-Keller reactor

- Temperature range: -45...250 °C
- Cooling power: 1.3 kW @ 250...0 °C
0.7 kW @ -20 °C
- Heating power: 1.5 kW / 3 kW
- Pump speed: 3300 rpm
- Hoses: 2x1 m; M24x1.5 (#9325)
- HTF: DW-Therm (#6479)
- Reactor: 1-litre jacketed glass reactor
- Reactor contents: 0.75 litre M90.055.03 (#6259)
- Reactor stirrer speed: 200 rpm
- Control: process

Unistat® 405w

Heating and cooling a Glas-Keller 1-litre glass reactor

Requirement

This case study examines the fast response of a Unistat 405w controlling the process temperature inside a 1-litre un-insulated glass reactor from the company "Glas-Keller" under two different control dynamics, "Fast-with overshoot" or "No overshoot".

Method

The Unistat 405w is connected to the Glas-Keller 1-litre reactor with two 1-metre insulated metal hoses. The reactor is filled with 0.75 litre of "M90.055.03", a silicon based HTF.

Results

The first cycle (20 °C to 60 °C to 20 °C) allows a small overshoot while the second cycle (20 °C to 60 °C to 20 °C) is "without overshoot" and so takes a longer time to reach set-point. It can be seen that even allowing for an overshoot, the control is so tight the overshoot in the first curve is negligible.

