



Setup details

Unistat® 425w & Buchi Glas Uster reactor

Temperature range: -40...250 °C

2.8 kW @ 250...100 °C Cooling power:

2.5 kW @ 0 °C 1.9 kW @ -20 °C 0.2 kW @ -40 °C

Heating power: 2.0 kW 2x1 m; M38x1.5 Hoses:

(#6656)

HTF: DW-Therm (#6479) Reactor: 20-litre jacketed glass

reactor

15 litre M90.055.03 Reactor content:

(#6259)

Stirrer: 150 rpm Control: process

Unistat® 425w

Heating a Buchi Glas Uster 20-litre jacketed glass reactor to 60 °C

Requirement

This case study looks at the performance of a Unistat 425w heating a 20-litre glass reactor from 20 °C to 60 °C under "process" control.

Method

The Unistat 425w is connected to the 20-litre Buchi Glas Uster reactor using two insulated metal 1-metre hoses. The reactor is filled with 15 litre of "M90.055.03", a silicon based HTF.

The jacket heats quickly to 102 °C creating a wide ΔT to pull the process temperature to its new set-point. As the process approaches the set-point the jacket cools rapidly to guide the process precisely to target temperature.

