

Unistat® 830

Cooling a Buchi Glas Uster 20-litre jacketed glass reactor to T_{min}

Requirement

This test is designed to determine the minimum achievable temperature that the Unistat 830 can cool a Buchi Glas Uster 20-litre reactor within 2 hours.

Method

The Unistat and reactor are connected using two 1.5-metre insulated metal hoses. The reactor is filled with 15 litre of "M90.055.03", a Huber supplied silicon based HTF.

Results

The "internal" (jacket) temperature cools to -57 °C in the first 20 minutes before the cooling power begins to asymptote. The final temperatures of jacket and process are -77 °C and -66 °C respectively.

Setup details

Unistat® 830 & Buchi Glas Uster reactor

Temperature range: -85...200 °C 3.6 kW @ 0 °C Cooling power: 2.2 kW @ -60 °C

3.6 @ 0 °C 3.5 @ -20...-40 °C 2.2 @ -60 °C 0.7 @ -80 °C

3 kW Heating power:

2x1.5 m; M38x1.5 Hoses:

(#6656)

HTF: DW-Therm (#6479) Reactor: 20-litre jacketed glass

reactor

15 litre M90.055.03 Reactor contents:

(#6259)

Reactor stirrer speed: 70 rpm Control: process



