



Setup details

Unistat® 830 & Buchi Glas Uster reactor

Temperature range: -85...200 °C 3.6 kW @ 0 °C Cooling power: 2.2 kW @ -60 °C

3.6 @ 0 °C 3.5 @ -20...-40 °C 2.2 @ -60 °C 0.7 @ -80 °C

Pump speed: 3500 rpm Heating power: 3 kW

2x1.5 m; M38x1.5 (#9616) Hoses: HTF: DW-Therm (#6479) Reactor: 20-litre jacketed glass

reactor

15 litre M90.055.03 Reactor contents:

(#6259)

Reactor stirrer speed: 70 rpm Control: process

Unistat® 830

Controlling simulated exothermic reactions in a Buchi Glas Uster 20-litre glass reactor at -40 °C

Requirement

This case study is to see the response of a Unistat 830 controlling simulated exothermic reactions in a Buchi Glas Uster 20-litre reactor at -40 °C:

- 1.) 50 Watt (43 kcal / hr)
- 2.) 100 Watt (86 kcal / hr)
- 3.) 150 Watt (129 kcal / hr)


Method

The Unistat and reactor are connected using two 1.5-metre insulated metal hoses. The reactor is filled with 15 litre of "M90.055.03", a Huber supplied silicon based HTF. The exothermic reactions are simulated using a controlled electric immersion heater.

Results

In each case the jacket temperature responds immediately to return the process temperature to its set-point and maintain exactly the setpoint during the "exotherm".


1.) 50 W (43 kcal / hr)







2.) 100 W (86 kcal / hr)



3.) 150 W (129 kcal / hr)

