

# Unistat® 510w

**Cooling and heating a Chemglass 50-litre jacketed glass reactor between 20 °C and -30 °C**

### Requirement

This case study demonstrates the speed of response and level of control achievable when a Chemglass 50-litre reactor is connected to a Unistat 510w.

### Method

The Unistat and reactor were connected using two 1.5 m insulated metal hoses. The reactor was filled with 37 litre of "M90.055.03", a Huber supplied silicon based HTF.

### Results

The process is ramped from 20 °C to -30 °C (50 K) in approximately 70 minutes while heating back to 20 °C takes approximately 25 minutes. In both cases the process reaches and maintains the new set-point with negligible over or under shoot.

### Setup details

Unistat® 510w & Chemglass 50-litre reactor

- Temperature range: -50...250 °C
- Cooling power: 5.3 kW @ 250...0 °C  
2.8 kW @ -20 °C  
0.9 kW @ -40 °C
- Heating power: 6.0 kW
- Hoses: 2x1.5 m; M30x1.5 (#6659)
- HTF: DW-Therm (#6479)
- Reactor: 50-litre Chemjacketed glass reactor (un-insulated)
- Reactor content: 37 litre M90.055.03
- Stirrer speed: 80 rpm
- Control: process

