

## Unistat® 510w

### Controlling simulated exothermic reactions within a Chemglass 50-litre glass reactor

#### Requirement

This case study shows the effectiveness of a Unistat 510w connected to a 50-litre reactor in the control of a three simulated exothermic reactions of differing strengths at 0 °C and 20 °C generated by an electric immersion heater.

#### Method

The Unistat and reactor were connected using two 1.5 m insulated metal hoses. The reactor was filled with 37 litre of "M90.055.03", a Huber supplied silicon based HTF.

#### Results

Once the "reaction" was under steady control the heater was turned "Off". It can be seen how rapidly the  $\Delta T$  between the process and jacket is increased to "suck" the thermal energy from the process to restore and maintain the process temperature set-point. The results can be viewed in the following graphics.


#### Setup details

Unistat® 510w & Chemglass 50-litre reactor


Temperature range:	-50...250 °C
Cooling power:	5.3 kW @ 250...0 °C 2.8 kW @ -20 °C 0.9 kW @ -40 °C
Heating power:	6.0 kW
Hoses:	2x1.5 m; M38x1.5 (#6659)
HTF:	DW-Therm (#6479)
Reactor:	50-litre Chemjacketed glass reactor (un-insulated)
Reactor content:	37 litre M90.055.03 (#6259)
Stirrer speed:	80 rpm
Control:	process



### 1. Simulated exothermic of 600 W at 0 °C



2. Simulated exothermic of 600 W at 20 °C



3. Simulated exothermic of 900 W at 20 °C

