

## Unistat® 910w

**Cooling a Diehm 100-litre jacketed glass reactor to -60°C**

### Requirement

This case study looks at the time taken for a Unistat 910w to cool the contents of a Diehm 100-litre un-insulated jacketed glass reactor to -60 °C from 20 °C.

### Method

The Unistat and reactor are connected using two 1.5-metre insulated metal hoses. The reactor is filled with 75 litre of "M90.055.03", a Huber supplied silicon based HTF.

A VPC-Bypass was installed to prevent a glass reactor breakage due to the pump pressure and reactor design.

### Results

It can be seen that though the 100-litre reactor represents a large thermal load to the Unistat 910w (the Unistat 910w is designed for use on reactors to a maximum of 50 litre), the process ramps through 80 K (20 °C to -60 °C) in around 140 minutes demonstrating great efficiencies in thermal transfer.

### Setup details

Unistat® 910w & Diehm 100-litre reactor

Temperature range: -90...250 °C  
 Cooling power: 5.2 kW from 250 °C to 20 °C  
 4.7 kW @ -40 °C  
 3.1 kW @ -60 °C  
 0.9 kW @ -80 °C  
 Heating power: 6.0 kW  
 Hoses: M38x1.5; 1x 2m #6657; 1x1m # 6655), VPC Bypass installed  
 HTF: M90.055.03 (#6259)  
 Reactor: 100-litre Diehm un-insulated jacketed glass reactor  
 Reactor content: 75 litre M90.055.03  
 Stirrer speed: 410 rpm  
 Control: process

