



Setup details

Stirrer speed:

Control:

Unistat® 510w & DDPS reactor

Temperature range: -50...250 °C

Cooling power: 5.3 kW @ 250...0 °C

2.8 kW @ -20 °C

0.9 kW @ -40 °C

Heating power: 6.0 kW

2x1.5 m; M38x1.5 Hoses:

(#6656)

HTF: DW-Therm (#6479) Reactor: 25-litre vacuum insu-

> lated jacketed glass reactor

Reactor content: 18.75 litre M90.055.03

(#6259) 80 rpm internal

Unistat® 510w

Heating a 25-litre DDPS reactor

Requirement

The graphic shows the performance of a Unistat 510w working with a 25-litre glass reactor connected together with M38x1.5 hoses.

Method

The DDPS 25-litre reactor and Unistat 425w were connected together with insulated "M38x1.5" hoses. The HTF circuit (reactor, Unistat and hoses) was filled with DW-Therm and the reactor was filled with 18.75 litre of "M90.055.03", a silicon based Huber supplied HTF.

The "internal" (jacket) temperature increases at a rate > 11 K per minute and reaches the set-point of 180 °C from -50 °C without any overshoot or undershoot within 21 minutes.

