

Ministat® 230-cc®-NR

Controlling a vacuum insulated Syrris 2-litre glass jacketed reactor to T_{min}

Requirement

This case study demonstrates the lowest achievable temperature, speed of cooling and heating and level of control when connected with a Syrris "Atlas" system configured with a 2-litre reactor.

Method

The reactor was filled to 1.4 litre with M90.055.03, the HTF used was Ethanol, the stirrer set to 700 rpm and the control to "process". The results were recorded using the "Spyware" software.

Results

It can be seen from the graphic that the Ministat 230-cc-NR cools the jacket to a minimum of -25 °C with a corresponding process temperature of approximately -24.5 °C within approximately 2 hours.

The heat up curve shows the precise control made possible by the Ministat 230-cc-NR as the process temperature reached exactly 20 °C from -25 °C in approximately 15 minutes.

Setup details

Ministat® 230-cc®-NR & Syrris vacuum insulated 2-litre glass jacketed reactor.

Temperature range: -40...200 °C Cooling power:

0.38 kW @ 0 °C 0.25 kW @ -20 °C

0.14 kW @ -30 °C

Heating power: 2 kW 4500 rpm Pump speed:

Hoses: 2x1 m; M16x1 (#9608)

HTF: Ethanol

2-litre jacketed glass Reactor:

reactor

1.4 litre M90.055.03 Reactor contents:

(#6259)

Reactor stirrer speed: 700 rpm

Control: process



