



CC[®]-505

CC[®]-505 cycling a 2 litre Radleys jacketed reactor

Requirement

This case study demonstrates the ability of the CC-505 refrigeration bath circulator to cycle the process temperature in a range from $+20^{\circ}$ C to -50° C, the closeness of the temperature control and the minimum process temperature achievable in the process mass.

Method

The 2 litre Radleys reactor was connected to the CC-505 using two M24x1.5 1-meter flexible hoses. The thermofluid used in the system was M40.165.10. "Process" control was carried out via a Pt100 sensor located in the process mass. Stirrer speed was set to 200 rpm.


Setup	details
Jetup	actunis

Temperature range:	-50°C+200°C
Cooling power:	1.2 kW @ +100°C
	1.2 kW @ +20°C
	1.0 kW @ 0°C
	0.6 kW @ -20°C
	0.15 kW @ -40°C
Heating power:	1.5 kW
Hoses:	M24x1.5; 2* 1 m
HTF:	M40.165.10
Reactor:	Radleys 2 litre
	jacketed reactor
Reactor content:	1.5 litre M40.165.10
Stirrer speed:	200 rpm
Control:	process

Results

Performance:

Cooling down and heating up in a range from +20°C to -20°C. The CC-505 needs approximately 45 minutes to cool down the reactor from +20°C to -20°C and 35 minutes to heat it up from -20°C to +20°C.



Lowest achievable temperature (T_{min}):

Once stable at +20°C under "Process" control, a set-point of -50°C is entered. The CC-505 cools the reactor down to the minimum achievable process temperature of -34°C.

