



## CC-510

**CC-510 cycling a 20-liter glass jacketed reactor**

### Requirement

This Case Study demonstrates the process temperature control and the minimum achievable process temperature when a CC-510 is connected with a Chemglass 20-liter glass jacketed reactor.

### Method

The 20-litre glass jacketed reactor was connected to CC-510 using M24 metal insulated hoses. The thermofluid used in the system was "M60.115/200.05". "Process" control was carried out via a Pt100 sensor located in the process mass. Stirrer speed was set to 100 rpm.


### Setup details

- Temperature range: -50°C...+200°C
- Cooling power: 2.1 kW @ +20°C
- 2.1 kW @ 0°C
- 1.0 kW @ -20°C
- Heating power: 3.0 kW
- Hoses: 2\*M24 metal insulated
- HTF: M60.115/200.05
- Reactor: Chemglass 20-liter glass jacketed reactor
- Reactor content: 15 l M60.115/200.05
- Stirrer speed: 100 rpm
- Control: process
- Amb. temperature: +25°C

## Results

### 1. Lowest achievable temperature (Tmin):

As the graphic shows, a minimum process temperature of -31.6°C was achieved.



## 2. Performance:

The table and graphic data show the speed and accuracy achieved as the CC-510 reaches each new set-point.

Start T	End T	Approximate time	Av. Ramp Rate
+20°C	-20°C	75 minutes	0.5 K/min
-20°C	+100°C	94 minutes	1.3 K/min
+100°C	+20°C	80 minutes	1.0 K/min

