



CC®-K6

CC®-K6 controlling 1 liter Asahi reactor

Requirement

This Case Study demonstrates the speed and accuracy when a CC-K6 is connected together with an Asahi 1 liter vacuum insulated reactor over the temperature range +20°C to +150°C and back to +20°C.

Method

The 1 liter Asahi vacuum insulated reactor, was connected to the CC-K6 using two metal insulated hoses. The thermofluid used in the system was M20.195/235.20. "Process" control was carried out via a Pt100 sensor located in the process mass. Stirrer speed was set to 150 rpm.

Setup details

Temperature range: -25°C...+200°C

Cooling power: 0,20 kW @ +20°C 0,15 kW @ 0°C

0,05 kW @ -20°C

 Heating power:
 2,0 kW

 Hoses:
 M16, 2 x 1 m

 HTF:
 M20.195/235.20

 Reactor:
 1 liter Asahi vacuum

insulated

Reactor content: 1 | M20.195/235.20


Stirrer speed: 150 rpm Control: Process Amb. temperature: +23°C

Results

1. Performance

The graphic shows the CC-K6 reaching and maintaining each new set point.

Start T	End T	Approximate time	Av. Ramp Rate	Fastest Ramp Rate
+20°C	+150°C	38 minutes	3.4 K/min	(+30°C to +60°C) 7.5 K/min
+150°C	+20°C	117 minutes	1.1 K/min	(+130°C to +100°C) 1.9 K/min





2. Lowest achievable temperature (Tmin):

The graphic shows the CC-K6 cooling the process to -7.8°C.

