

Unistat® 705w

Heating and cooling a DDPS 2-litre jacketed reactor

Requirement

The purpose of this case study is to demonstrate the performance of a Unistat 705w in heating and cooling a 2-litre glass reactor.

Method

The Unistat and reactor are connected using two 1-metre insulated metal hoses. The reactor is filled with 1.5 litre of "M90.055.03", a Huber supplied silicon based HTF.

Result

The resulting curve shows it takes 33 minutes to heat from 20 °C to 100 °C giving a heating ramp-rate of 2.4 K/min.

In cooling, the curve shows the process takes 35 minutes to reach 20 °C giving a cooling ramp-rate of 2.3 K/min.

Setup details

Heating power:

Unistat® 705w & 2-litre DDPS reactor

Temperature range: -75...250 °C

0.6 kW @ 250...100 °C Cooling power:

0.65 kW @ 0 °C 0.6 kW @ -20...-40 °C 0.3 kW @ -60 °C

1.5 kW/3 kW (in this test

limited to 1.0 kW)

Pump speed: 3300 rpm

2x1 m; M24x1.5 (#9325) Hoses: HTF: DW-Therm (#6479) Reactor: 2-litre un-insulated

iacketed glass reactor

Reactor content: 1.5 litre M90.055.03

(#6259)

Stirrer speed: 200 rpm Control: process



