



Unistat[®] 910w

Heating a 50-litre jacketed glass reactor from -50 °C to 20 °C

Requirement

This case study demonstrates the response of Unistat 910w to heat the contents of an uninsulated 50-litre glass reactor from -50 °C to 20 °C.

Method

The Unistat and reactor are connected using two 1.5-metre insulated metal hoses. The reactor is filled with 35 litre of "M90.055.03", a Huber supplied silicon based HTF.

Results

The process temperature is ramped rapidly to set-point with a negligible overshoot and achieving stability within 40 minutes.

Setup details

Unistat[®] 910w & Chemglass reactor

Temperature range:	-90250 °C
Cooling power:	5.2 kW @ 25020 °C
	4.7 kW @ -40 °C
Heating power:	6.0 kW
Hoses:	2x1.5 m; M30x1.5
	(#6386)
HTF:	DW-Therm (#6479)
Reactor:	50-litre un-insulated
	jacketed glass reactor
Reactor content	35 litre M90.055.03
	(#6259)
Stirrer speed	80 rpm
Control	process



