

Unistat® 610w

Control dynamics on a Buchi Glas Uster 20-litre reactor

Requirement

Every Unistat can be set to ramp "Fast with small overshoot" or "No overshoot". This case study looks at the response under different "control dynamics" of a Unistat 610w when cooling and heating a Buchi Glas Uster 20-litre reactor between 20 °C and -60 °C.

Method

The Unistat and reactor were connected using two 1.5 m insulated metal hoses. The reactor was filled with 15 litre of "M90.055.03", a Huber supplied silicon based HTF.

Results

The first and second curves (20 °C to 60 °C and back to 20 °C) show the function of "Periodic - fast, small overshoot" control dynamic. The first curve illustrates a heating process in a time of 16 minutes. It can be clearly seen that the internal temperature heats to 103 °C, thus the process temperature reaches -60 °C very quickly. The Unistat

610w cools the 20-litre reactor back to 20 °C in approximately 17 minutes through a ΔT of 40 K.

The third curve (20 °C to 60 °C) shows the same temperature profile but with "Aperiodic - no overshoot" control. The Unistat takes slightly longer (30 minutes) to heat to avoid any over or undershoot of the set-point.

Setup details

Unistat® 610w & Buchi Glas Uster reactor

| | |
|--------------------|--|
| Temperature range: | -60...200 °C |
| Cooling power: | 7.0 kW @ 200...0 °C 6.4 kW @ -20 °C 3.3 kW @ -40 °C 0.8 kW @ -60 °C |
| Heating power: | 6.0 kW |
| Hoses: | 2x1.5 m; M38x1.5 (#6656) |
| HTF: | DW-Therm (#6479) |
| Reactor: | 20-litre jacketed glass reactor |
| Reactor content: | 15 litre DW-Therm (#6479) |
| Stirrer speed: | 70 rpm |
| Control: | process |

