



Unistat[®] 825w

Heating and cooling a Buchi Glas Uster 10-litre jacketed reactor

Requirement

This case study intended to investigate the performance of a Unistat 825w heating and cooling a Buchi Glas Uster 10-litre reactor between 20 °C to 100 °C and then from 20 °C to 180 °C.

Method

The Unistat and reactor are connected using two 1.5-metre insulated metal hoses. The reactor is filled with 7.5 litre of "M90.055.03", a Huber supplied silicon based HTF.

Results

In the first segment (20 °C to 100 °C) the heating ramp rate of 2.94 K/min. heats the process temperature to 100 °C in 47 minutes. In the second segment (20 °C to 180 °C) the average heating ramp rate of 1.9 K/min. brings the process temperature to the set-point within 1:47 hour.

Setup details

Unistat® 825w & Buchi «miniPilot» 10 reactor

Temperature range:	-85250 °C
Cooling power:	2.4 kW @ 040 °C
	1.5 kW @ -60 °C
Heating power:	3.0 kW
Pump speed:	3500 rpm
Hoses:	2x1.5 m; M30x1.5
	(#6386)
HTF:	DW-Therm (#6479)
Reactor:	10-litre jacketed glass
	reactor
Reactor contents:	7.5 litre M90.055.03
	(#6259)
Reactor stirrer speed:	400 rpm
Control:	process



